




MICRO-TECH™
ENDOSCOPY


TRIPLE-LUMEN SPHINCTEROTOME
PAPILLOTOMES

Tapered Tip
Monofilament cutting wire



Your Trusted Partner In Endoscopy!

TRIPLE-LUMEN SPHINCTEROTOME



SPECIFICATIONS

REF	Lumen	Tip Length (mm)	Cut Wire Length (mm)	Working length (mm)	Recommended Guidewire (inch)	Partly Insulated Cut Wire
DSP-30220-111111	3	2	20	2000	0.035	NO
DSP-30220-111211	3	2	20	2000	0.035	YES
DSP-30225-111111	3	2	25	2000	0.035	NO
DSP-30225-111211	3	2	25	2000	0.035	YES
DSP-30230-111111	3	2	30	2000	0.035	NO
DSP-30230-111211	3	2	30	2000	0.035	YES
DSP-30520-111111	3	5	20	2000	0.035	NO
DSP-30520-111211	3	5	20	2000	0.035	YES
DSP-30525-111111	3	5	25	2000	0.035	NO
DSP-30525-111211	3	5	25	2000	0.035	YES
DSP-30530-111111	3	5	30	2000	0.035	NO
DSP-30530-111211	3	5	30	2000	0.035	YES

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National
Hi-Tech Zone, Nanjing, China, 210032
+86-25-58646393

Micro-tech.com.cn/en/

Sales@micro-tech.com.cn