



MICRO-TECH™

ENDOSCOPY

SINGLE-USE ELECTROSURGICAL

## ESD KNIFE

Provides for marking, elevation, incision/dissection and coagulation in a single instrument.

HIGH-PRESSURE WATERJET



YOUR TRUSTED PARTNER IN ENDOSCOPY!

## KEY BENEFITS

- Integrated electrosurgical and waterjet functions, eliminating the need for instrument change and saving valuable procedure time.
- High-pressure waterjet creates submucosal water cushion to prevent thermal and mechanical injury during resection, ensuring maximum safety.
- Three types of electrosurgical tip available for various clinical applications.

I-type: Versatile application.

T-type: Marking/dissection and less likely to slip.

O-type: Spherical tip maximizes coagulation effect.



I-type



T-type



O-type

## SPECIFICATIONS

REF	Shape of the Cutting Knife	Cutting Knife Length (mm)	Working Length (mm)	Channel Size Ø (mm)
MK-T-1-195	T	1.5	1950	≥2.8
MK-T-2-195	T	2	1950	≥2.8
MK-T-4-195	T	4	1950	≥2.8
MK-I-1-195	I	1.5	1950	≥2.8
MK-I-2-195	I	2	1950	≥2.8
MK-I-4-195	I	4	1950	≥2.8
MK-O-1-195	O	1.5	1950	≥2.8
MK-O-2-195	O	2	1950	≥2.8
MK-O-4-195	O	4	1950	≥2.8

1650mm and 2350mm working lengths optional

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.