



LOCKADO™ CLIP

NEW

STRONGER HOLD, MORE COMFORT

The new LOCKADO™ Clip from MICRO-TECH is designed to control bleeding more efficiently than ever. This is achieved via ten teeth that operate as anchoring spikes to ensure optimal fixation on the tissue of the bleeding area. Users can imitate the clip's actual closure while being able to reopen it at will. So it's repositionable and this ensures more precision and flexibility which in turn reduces the total operational time. We have also increased the compression pressure by 15%, which

enables a significantly improved closing force. For maximum patient comfort, the clip segment projecting into the lumen was designed to be particularly short. This enables comfortable nutrition and minimizes the risk of injury to the opposing tract. MICRO-TECH offers the LOCKADO™ in different options. For example you can choose between three clip widths (8, 11 and 16 mm) to cover all needs and possible applications.

SPECIFIC CHARACTERISTICS

- Bleeding is better targeted and can be stopped quickly and specifically
- Ten teeth for ideal hold
- Clip arms with maximum contact pressure
- Unlimited repositioning possible prior to fixation
- High fixation due to better adaptation of tissue
- Increased the clips' compression pressure by 15%
- Short clip segment for greater patient comfort



Maximum span width of 16 mm

SPECIFICATIONS

REF	Clip max. span width mm	Length mm	Instrument Ø mm	Grip color	Covering
LOCKADO™					
LOCK-C-26-195-C	8	1950	2.6	Pink	Pink
LOCK-C-26-230-C	8	2300	2.6	Pink	Pink
LOCK-D-26-195-C	11	1950	2.6	Blue	Blue
LOCK-D-26-230-C	11	2300	2.6	Blue	Blue
LOCK-F-26-195-C	16	1950	2.6	Yellow	Yellow
LOCK-F-26-230-C	16	2300	2.6	Yellow	Yellow
LOCKADO™ PREMIUM-ECO					
LOCK-C-26-195-C-S	8	1950	2.6	Pink	Grey
LOCK-C-26-230-C-S	8	2300	2.6	Pink	Grey
LOCK-D-26-195-C-S	11	1950	2.6	Blue	Grey
LOCK-D-26-230-C-S	11	2300	2.6	Blue	Grey
LOCK-F-26-195-C-S	16	1950	2.6	Yellow	Grey
LOCK-F-26-230-C-S	16	2300	2.6	Yellow	Grey
LOCKADO™ ECO					
LOCK-C-26-195	8	1950	2.6	Pink	Without
LOCK-C-26-230	8	2300	2.6	Pink	Without
LOCK-D-26-195	11	1950	2.6	Blue	Without
LOCK-D-26-230	11	2300	2.6	Blue	Without
LOCK-F-26-195	16	1950	2.6	Yellow	Without
LOCK-F-26-230	16	2300	2.6	Yellow	Without

Packaging unit: 10 pieces

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need