




MICRO-TECH™  
ENDOSCOPY

SINGLE-USE COLD SNARE  
**DIAMOND CUT™**  
Speed and precision with every cut

THIN CUTTING WIRE




10 OR 15MM DIAMOND LOOP

Your **Trusted Partner** In Endoscopy!

## KEY BENEFITS

- Diamond-shaped wire is specifically designed for resecting diminutive polyps while using a cold snaring technique.
- The thin 0.22mm wire provides a “guillotine-type” cut.
- Stiffer catheter aids in insertion through an extended scope position and provides support during resection.
- Cold snaring offers several clinical advantages for polyps under 1cm in diameter.
- Comfortable ergonomic handle accommodates various hand sizes.



## SPECIFICATIONS

REF	Loop Shape	Loop Diameter (mm)	Working Length (cm)	Min. Working Channel(mm)
CS3-11023230	Diamond	10	230	2.8
CS3-11523230	Diamond	15	230	2.8

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National  
Hi-Tech Zone, Nanjing, China, 210032  
+86-25-58646393

[Micro-tech.com.cn/en/](http://Micro-tech.com.cn/en/)

[Sales@micro-tech.com.cn](mailto:Sales@micro-tech.com.cn)