




## FREEDOM FROM REPROCESSING **DISPOSABLE ACCESSORIES FOR ESOPHAGEAL DILATORS**

Repeated use and reprocessing of reusable guidewires and cleaning brushes can result in worn, frayed, bent or otherwise damaged devices.


The use of disposable guidewires and disposable cleaning brushes can reduce the risk of device damage as well as the potential cross-contamination from using reprocessed accessories with reusable polyvinyl dilator systems.

The Single-Use Marked Spring Tip Guidewire and Dilator Cleaning Brush are for use with American and Savary dilators.




**YOUR TRUSTED PARTNER IN ENDOSCOPY!**

## SPECIFICATIONS

| REF                                | Diameter Of Spring Soft (mm) | Working Length (mm) | Diameter Of Guide Wire Rod (mm) | Packaging   |
|------------------------------------|------------------------------|---------------------|---------------------------------|---|
| <b>Marked Spring Tip Guidewire</b> |                              |                     |                                 |  |
| DG-18-08-21                        | 1.8                          | 2100                | 0.8                             | 10/Box  |



| REF                           | Diameter Of Catheter (mm) | Working Length (mm) | Diameter Of Brush Head (mm) | Packaging   |
|-------------------------------|---------------------------|---------------------|-----------------------------|---|
| <b>Dilator Cleaning Brush</b> |                           |                     |                             |  |
| CB06-0040/09-S                | 1                         | 900                 | 3                           | 10/Box  |
| CB10-0030/09-S                | 0.6                       | 900                 | 4                           | 10/Box  |

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.