



MICRO-TECHTM
ENDOSCOPY

SINGLE-USE HEMOCLIP FAMILY

ECO CLIP

Repeated opening and closing possible
Wide opening
Excellent holding power
Smooth running under tortuous anatomy



YOUR TRUSTED PARTNER IN ENDOSCOPY!

CLIP WITH CONFIDENCE

The innovative clipping system from MICRO-TECH, combines precise positioning with a reliable holding force. This makes the system suitable for effective hemostasis as well as for endoscopic markings. The clip can be opened and closed multiple times, thus facilitating the exact repositioning of the clip. The clip can be used with any endoscopes equipped with working channels with a minimum diameter of 2.8mm.



SPECIFICATIONS

REF	Clip max. span width mm	Length mm	Instrument Ø (mm)	Grip color	Covering
ROCC-C-26-195	8	1950	2.6	Pink	Without
ROCC-C-26-230	8	2300	2.6	Pink	Without
ROCC-D-26-195	11	1950	2.6	Blue	Without
ROCC-D-26-230	11	2300	2.6	Blue	Without
ROCC-F-26-195	16	1950	2.6	Yellow	Without
ROCC-F-26-230	16	2300	2.6	Yellow	Without

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need