

Injection Needles



Performance Characteristics:

- Lock to ensure that it would not accidentally extend during insertion to damage the scope and to have a stable needle penetration.
- Needle release button design, gently press the button to return the needle fully into the sheath for protecting the endoscope from any damage.
- Ergonomic handle design: complete needle insertion and with drawal with one hand.



REF	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
IN02-224231802	22G	4	2.3	180	10/Box
IN02-224232002	22G	4	2.3	200	10/Box
IN02-224232302	22G	4	2.3	230	10/Box
IN02-225231802	22G	5	2.3	180	10/Box
IN02-225232002	22G	5	2.3	200	10/Box
IN02-225232302	22G	5	2.3	230	10/Box
IN02-226231802	22G	6	2.3	180	10/Box
IN02-226232002	22G	6	2.3	200	10/Box
IN02-226232302	22G	6	2.3	230	10/Box
IN02-254231802	25G	4	2.3	180	10/Box
IN02-254232002	25G	4	2.3	200	10/Box
IN02-254232302	25G	4	2.3	230	10/Box
IN02-255231802	25G	5	2.3	180	10/Box
IN02-255232002	25G	5	2.3	200	10/Box
IN02-255232302	25G	5	2.3	230	10/Box
IN02-256231802	25G	6	2.3	180	10/Box
IN02-256232002	25G	6	2.3	200	10/Box
IN02-256232302	25G	6	2.3	230	10/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National
Hi-Tech Zone, Nanjing, China, 210032
+86-25-58646393

Micro-tech.com.cn/en/

Sales@micro-tech.com.cn