




MICRO-TECHTM
ENDOSCOPY

FULLY COVERED ESOPHAGEAL STENTS


SEGMENTED STENT



Your **Trusted Partner** In Endoscopy!

FLEXIBILITY REDEFINED

- Segmented modular stent body allowing independent movements between each module to maximize conformance to curved or tortuous anatomy.
- Segmented design minimizes migration
- Comparable radial force
- Corrosion resistant and elastic covering
- Retrieval loop for removal and repositioning



Esophageal Stents with segmented design to allow higher flexibility and conformity for bypassing esophageal stenoses.

SPECIFICATIONS

Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Stent Length (mm)	Delivery System (mm)	Working length (mm)
16	22	60/80/100/120/140	6.7	650
18	24			
20	26			
22	28			

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National
Hi-Tech Zone, Nanjing, China, 210032
+86-25-58646393

Micro-tech.com.cn/en/

Sales@micro-tech.com.cn