



# SHORT WIRE METAL STENT

## FOR SAFE AND EASY STENT RELEASE

- Compatible with our short wire guide wires
- Platinum Line for excellent visibility
- Safe application
- Easier release of the stent without the need for an irrigation port
- Cylindrical stentform

# SPECIFICATIONS

REF	Ø mm	Length mm	Covering mm	Introducer max Ø F	Working length mm	Resheathable	Recommended guidewire inch
<b>SHORT WIRE METAL STENTS</b>							
RST43-001-10.040	10	40	0	8.5	1900	Yes	0.035"
RST43-001-10.060	10	60	0	8.5	1900	Yes	0.035"
RST43-001-10.080	10	80	0	8.5	1900	Yes	0.035"
RST43-001-10.100	10	100	0	8.5	1900	Yes	0.035"
RST43-002-10.040	10	40	30	8.5	1900	Yes	0.035"
RST43-002-10.060	10	60	50	8.5	1900	Yes	0.035"
RST43-002-10.080	10	80	70	8.5	1900	Yes	0.035"
RST43-002-10.100	10	100	90	8.5	1900	Yes	0.035"
RST43-004-10.040	10	40	40	8.5	1900	No	0.035"
RST43-004-10.060	10	60	60	8.5	1900	No	0.035"
RST43-004-10.080	10	80	80	8.5	1900	No	0.035"
RST43-004-10.100	10	100	100	8.5	1900	No	0.035"

REF	Ø inches	Total length mm	Tip length mm	Tip form	Set-up	Covering
<b>ECO GUIDE WIRE</b>						
MTN-BM-89/26-A	0.035"	2600	65	straight	hydrophilic, radiopaque	blue-yellow covering
MTN-BM-89/26-A-J	0.035"	2600	65	J-tip	hydrophilic, radiopaque	blue-yellow covering
MTN-BM-63/26-A	0.025"	2600	65	straight	hydrophilic, radiopaque	blue-yellow covering
MTN-BM-63/26-A-J	0.025"	2600	65	J-tip	hydrophilic, radiopaque	blue-yellow covering

Packaging unit Eco Guide Wire: 2 pieces

subject to errors and technical alterations. Rev 11/2019