



**MICRO-TECH**<sup>TM</sup>  
ENDOSCOPY

# POLYPECTOMY SNARES


Speed and precision with every cut



**Your Trusted Partner In Endoscopy!**

## KEY BENEFITS

- Very comfortable to use
- Reliable and safe
- Innovative power connection
- Ergonomic 3-ring handle
- Excellent cutting characteristics
- High strength of shape



## SPECIFICATIONS

### • Snares - Stiff

REF	Snare Shape	Snare Diameter (mm)	Sheath O.D. (mm)	Channel Size (mm)	Working Length (cm)	Packaging
PFS01-01023230	Oval	10	2.3	≥φ2.8	230	10/Box
PFS01-01523230	Oval	15	2.3	≥φ2.8	230	10/Box
PFS01-02023230	Oval	20	2.3	≥φ2.8	230	10/Box
PFS01-02423230	Oval	24	2.3	≥φ2.8	230	10/Box
PFS01-03023230	Oval	30	2.3	≥φ2.8	230	10/Box
PFS01-03623230	Oval	36	2.3	≥φ2.8	230	10/Box
PFS01-11523230	Hexagonal	15	2.3	≥φ2.8	230	10/Box
PFS01-12423230	Hexagonal	24	2.3	≥φ2.8	230	10/Box
PFS01-33023230	Combinational	30	2.3	≥φ2.8	230	10/Box
PFS01-01023230	Oval	10	2.3	≥φ2.8	230	10/Box

### • Snares - Flexible

REF	Snare Shape	Snare Diameter (mm)	Sheath O.D. (mm)	Channel Size (mm)	Working Length (cm)	Packaging
PFS02-01023230	Oval	10	2.3	≥φ2.8	230	10/Box
PFS02-01523230	Oval	15	2.3	≥φ2.8	230	10/Box
PFS02-02023230	Oval	20	2.3	≥φ2.8	230	10/Box
PFS02-02423230	Oval	24	2.3	≥φ2.8	230	10/Box
PFS02-03023230	Oval	30	2.3	≥φ2.8	230	10/Box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National  
Hi-Tech Zone, Nanjing, China, 210032  
+86-25-58646393

[Micro-tech.com.cn/en/](http://Micro-tech.com.cn/en/)

[Sales@micro-tech.com.cn](mailto:Sales@micro-tech.com.cn)