




**MICRO-TECH**<sup>TM</sup>  
**ENDOSCOPY**


Multi-Stage Stone Extraction Balloon

# **SWEEPER**<sup>TM</sup>

Inflated balloon presents cylindrical shape  
Tapered Tip



High-quality latex  
3-stage balloons






3 pre-measured syringes  
for reliable inflation

Stiffer catheter with  
excellent pushability

**Your Trusted Partner In Endoscopy!**

# MULTI-STAGE STONE EXTRACTION BALLOON



## SPECIFICATIONS

REF	Balloon Inflated O.D (mm)	Working Length (cm)	Recommended Guidewire (in)	Minimum Channel Size (mm)
SRB-T9/12/15-20	9-12-15	200	0.035	2.8
SRB-T12/15/18-20	12-15-18	200	0.035	3.2
SRB-T15/18/20-21	15-18-21	200	0.035	3.2

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National  
Hi-Tech Zone, Nanjing, China, 210032  
+86-25-58646393

[Micro-tech.com.cn/en/](http://Micro-tech.com.cn/en/)

[Sales@micro-tech.com.cn](mailto:Sales@micro-tech.com.cn)