



**MICRO-TECH**<sup>TM</sup>  
**ENDOSCOPY**

SINGLE-USE HEMOCLIP FAMILY

# **SURECLIP**<sup>TM</sup>

Repositioning  
Reliable 360° Rotation  
Shorter Stem



**YOUR TRUSTED PARTNER IN ENDOSCOPY**

# CLIP WITH CONFIDENCE

## REPOSITIONING

SureClip's unique design permits repositioning prior to deployment. Being able to reposition a clip, may help improve placement accuracy.

## RELIABLE 360° ROTATION

SureClip can be rotated, which helps to provide the correct orientation for tissue approximation or vessel closure.

## SHORTER STEM

A shorter stem makes the clip less obtrusive, improving visualization of the target area, particularly when multiple clips are placed in close proximity.



## SPECIFICATIONS

	REF	Clip Opening(mm)	Working Length(cm)	Min Working Channel(mm)
SureClip™ mini	ROCC-C-26-165-C	8	165	2.8
	ROCC-C-26-195-C	8	195	2.8
	ROCC-C-26-230-C	8	230	2.8
SureClip™	ROCC-D-26-165-C	11	165	2.8
	ROCC-D-26-195-C	11	195	2.8
	ROCC-D-26-230-C	11	230	2.8
SureClip™ plus	ROCC-F-26-165-C	16	165	2.8
	ROCC-F-26-195-C	16	195	2.8
	ROCC-F-26-230-C	16	230	2.8

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.