



MICRO-TECHTM
ENDOSCOPY

Multi-Stage Stone Extraction Balloon

SWEEPERTM

Inflated balloon presents cylindrical shape
Tapered Tip



High-quality latex
3-stage balloons

3 pre-measured syringes
for reliable inflation

Stiffer catheter with
excellent pushability

Your **Trusted Partner** In Endoscopy!

MULTI-STAGE STONE EXTRACTION BALLOON



SPECIFICATIONS

REF	Balloon Inflated O.D (mm)	Working Length (cm)	Recommended Guidewire (in)	Minimum Channel Size (mm)
SRB-T9/12/15-20	9-12-15	200	0.035	2.8
SRB-T12/15/18-20	12-15-18	200	0.035	3.2
SRB-T15/18/20-21	15-18-21	200	0.035	3.2

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



No.10, Gaoke 3rd Road, National
Hi-Tech Zone, Nanjing, China, 210032
+86-25-58646393

Micro-tech.com.cn/en/

Sales@micro-tech.com.cn