



TECHBITE™ BIOPSY FORCEPS

FOR CLEAN AND LARGE BIOPSATES

For large and precisely cut biopsies, the Techbite™ forceps from MICRO-TECH are equipped with a special design. The spoons allow a meaningful sample capacity and are characterized by a precise cutting performance. To meet the different clinical requirements the assort-

ment includes forceps with and without spike, as well as with smooth and serrated edge. The thumb ring handle offers optimal control and a clear tactile feeling. In order to protect the working channel and reduce sliding friction, the Techbite™ is completely covered.

SPECIFIC CHARACTERISTICS

- Large spoon design and precise cuts
- With smooth and serrated cutting edge
- Available with or without spike
- Complete coating reduces sliding friction
- Thumb ring handle for optimal control



SPECIFICATIONS

REF	Jaw	Span width jaw mm	Jaw Ø mm	Spike	Working length mm	Covering	Working channel Ø mm
TECHBITE™ BIOPSY FORCEPS							
EBF33-11018120	standard	4.5	1.8	No	1200	blue	2.0
EBF33-11018180	standard	4.5	1.8	No	1800	blue	2.0
EBF33-11018230	standard	4.5	1.8	No	2300	blue	2.0
EBF33-11023120	standard	6.7	2.3	No	1200	blue	2.8
EBF33-11023180	standard	6.7	2.3	No	1800	blue	2.8
EBF33-11023230	standard	6.7	2.3	No	2300	blue	2.8
EBF33-11123120	standard	6.7	2.3	Yes	1200	blue	2.8
EBF33-11123180	standard	6.7	2.3	Yes	1800	blue	2.8
EBF33-11123230	standard	6.7	2.3	Yes	2300	blue	2.8
EBF43-11018120	serrated	4.5	1.8	No	1200	blue	2.0
EBF43-11018180	serrated	4.5	1.8	No	1800	blue	2.0
EBF43-11018230	serrated	4.5	1.8	No	2300	blue	2.0
EBF43-11023120	serrated	6.7	2.3	No	1200	blue	2.8
EBF43-11023180	serrated	6.7	2.3	No	1800	blue	2.8
EBF43-11023230	serrated	6.7	2.3	No	2300	blue	2.8
EBF43-11123120	serrated	6.7	2.3	Yes	1200	blue	2.8
EBF43-11123180	serrated	6.7	2.3	Yes	1800	blue	2.8
EBF43-11123230	serrated	6.7	2.3	Yes	2300	blue	2.8

Packaging unit: 10 pieces

Subject to errors and technical alterations. Rev.09/2019